

Procedure Pack

UroFunnel® Drape pack with Trolley Cover

Description

- Article 12 Sterile Procedure Pack for Urology
- Tyvek pouch/Headerbag
- Wrapped in Trolley Cover
- works with any UroFunnel Frame
- 3 years life
- No IFU Required
- Latex Free

Product Ref: **PM-UP-8002TC**



Intended Use: For intended Use see procedure pack component specifications

Ergonomic considerations: N/A

Storage and Handling: Store in a cool dry environment away from direct sunlight at room temperature between 10 and 25 C.

Pack Contents:

QTY	REF	Description	CE
1	TC-01	Trolley cover 150x200cm	CE
1	PM-1052	UroFunnel Bag	CE
1	PM-5100	UroFunnel Patient Drape	CE

Pack Specification: PM-8002-TC, Issue 19-2, 20190930