

Procedure Pack

Product Ref: **PM-UP-8003**

UroFunnel® Cystoscopy pack

Description

- Article 12 Sterile Procedure Pack for Urology
- Tyvek pouch/Headerbag
- Wrapped in Trolley Cover
- works with any UroFunnel Frame
- 3 years life
- No IFU Required
- Latex Free



Intended Use: For intended Use see procedure pack component specifications

Ergonomic considerations: N/A

Storage and Handling: Store in a cool dry environment away from direct sunlight at room temperature between 10 and 25 C.

Pack Contents:

QTY	REF	Description	CE
1	TC-01	Trolley cover 150x200cm	CE
1	PM-1052	UroFunnel Bag	CE
1	PM-5100	UroFunnel Patient Drape	CE
1	KD-800	Kidney Dish, 800ml	CE
1	B-120	Galipot, 120ml	CE
1	B-250-B	Galipot, 250ml	CE
1	SYR-20LS	Syringe, Luer Slip, 20ml	CE0050
1	9002	Camera Drape, Ring, Elastic	CE
1	SW-10XR	Swabs, Gauze, 10x10cm X-Ray	CE0123
2	SH-18	Sponge Holders, 18.75cm	CE
1	9502	Irrigation Set, Double, Cuff and Spike	CE2797

Pack Specification: PM-8003, Issue 19-2, 20190930