

## Procedure Pack

Product Ref: **PM-UP-8526 PCNL**

PCNL pack

### Description

- Article 12 Sterile Procedure Pack for Urology
- Tyvek pouch/Headerbag
- Wrapped in Trolley Cover
- works with any UroFunnel Frame
- 3 years life
- No IFU Required
- Latex Free



**Intended Use:** For intended Use see procedure pack component specifications

**Ergonomic considerations:** N/A

**Storage and Handling:** Store in a cool dry environment away from direct sunlight at room temperature between 10 and 25 C.

### Pack Contents:

QTY	REF	Description	CE
1	TC-01	Trolley cover 150x200cm	CE
2	PM-1052	UroFunnel Bag	CE
1	PM-5301	Nephrostomy Drape	CE
1	AT-01	Absorbent Towels	CE
1	HT-01	Hand Towels	CE
1	TRAY-01	Tray, 33x23x5.5, Blue	CE
1	GW-10	Guide Wire Bowl	CE
1	KD-800	Kidney Dish, 800ml	CE
1	B-120	Galipot, 120ml	CE
1	B-250-B	Galipot, 250ml, Blue	CE
1	SYR-20LS	Syringe, Luer Slip, 20ml	CE0050
1	9003	Camera Drape, Closed, Spring, Sleeve	CE
1	9000R	Camera Drape, Tape, Ring	CE
1	SW-10XR	Swabs, Gauze, 10x10cm X-Ray	CE0123
2	BD-32	Bull dog clip, 32mm	<i>Not a medical device</i>
2	SH-18	Sponge Holders, 18.75cm	CE
1	9502	Irrigation Set, Double, Cuff and Spike	CE2797
1	SC-13S	Scissors, 13cm, Sharp	CE
1	FC-11FT	Forceps, Fine Tip	CE

Pack Specification: PM-8526, Issue 19-2, 20190930